

Overview

The T-Cup panel drug test detects the most commonly abused drugs. The OneStep cup is manufactured for accuracy. Includes "O" ring in the cap to prevent leakage. The T-Cup All-In-One Drug Test Cup can be used for the synchronal and qualitative detection of most of the drug metabolites found in human urine at particular cutoff levels. It supplies quick results with an easy-to-read visual panel.

TEST	CODE
Amphetamine	AMP
Barbiturates	BAR
Benzodiazepines	BZO
Buprenorphine	BUP
Cocaine	COC
Synthetic Cannabinoids	K2/SPICE
Marijuana	THC
Methadone	MTD
Methamphetamine	M-AMP
Opiates	OPI
Alcohol	ETG
Phencyclidine	ETG
Tricyclic Antidepressants	TCA
Ecstasy	MDMA
Oxycodone	OXY
Propoxyphene	PPX
Fentanyl	FTY
Tramadol	TRAM

Minimum 12-month shelf life.

Attributes

Option	AMP, BUP, BZO, COC, M-AMP, MDMA, MTD, OPI, K2/SPICE, ETG, FTY, TRAM, OXY, THC, AMP, BAR, BUP, BZO, COC, ETG, M-AMP, MDMA, MTD, OPI, K2/SPICE, ETG, FTY, TRAM, OXY, THC
Packed	25/Cs
Product Unit Type	Case

Images

